

Product Datasheet

SLC12A2 antibody (orb631222)

Description	Rabbit polyclonal antibody to SLC12A2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	solute carrier family 12(sodium/potassium/chloride transporters), member 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	66 kDa
Uniprot ID	P55011
Dilution Range	WB: 1:500-1:5000, IP: 1:200-1:2000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.