



Product Datasheet

SNAP25 antibody (orb631358)


Description	Mouse monoclonal antibody to SNAP25
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	synaptosomal-associated protein, 25kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	28kd
Uniprot ID	P60880
Dilution Range	WB: 1:5000-1:20000, IP: 1:500-1:1000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of mouse brain tissue



Immunofluorescent analysis of PC-12 cells



Western blot analysis of HEK-293 cells