

Product Datasheet

SOD2 antibody (orb631400)

Description	Rabbit polyclonal antibody to SOD2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	superoxide dismutase 2, mitochondrial
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	16 kDa
Uniprot ID	P04179
Dilution Range	WB: 1:500-1:5000, IP: 1:500-1:5000, IHC-P: 1:100-1:400, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.