

Product Datasheet

SRC antibody (orb631497)

Description	Rabbit polyclonal antibody to SRC
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	v-src sarcoma(Schmidt-Ruppin A-2) viral oncogene homolog(avian)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	55kd,60kd
Uniprot ID	P12931
Dilution Range	WB: 1:500-1:2000, IP: 1:500-1:2000, IHC-P: 1:100-1:400, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.