

Product Datasheet

SRP14 antibody (orb631506)

Description	Rabbit polyclonal antibody to SRP14
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	signal recognition particle 14kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	18kd
Uniprot ID	P37108
Dilution Range	WB: 1:500-1:1000, IP: 1:200-1:1000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.