

Product Datasheet

TAF3 antibody (orb631709)

Description	Rabbit polyclonal antibody to TAF3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	TAF3 RNA polymerase II, TATA box binding protein(TBP)-associated factor, 140kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	80 kDa,130 kDa
Uniprot ID	Q5VWG9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.