

Product Datasheet

TCIRG1 antibody (orb631779)

Description	Rabbit polyclonal antibody to TCIRG1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	T-cell, immune regulator 1, ATPase, H+ transporting, lysosomal V0 subunit A3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	65-70 kDa
Uniprot ID	Q13488
Dilution Range	WB: 1:500-1:2000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.