

Product Datasheet

TEAD1 antibody (orb631797)

Description	Rabbit polyclonal antibody to TEAD1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	TEA domain family member 1 (SV40 transcriptional enhancer factor)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	P28347
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, CHIP: 1:20-1:50
Expiration Date	12 months from date of receipt.