

Product Datasheet

TFAM antibody (orb631828)

Description	Rabbit polyclonal antibody to TFAM
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	transcription factor A, mitochondrial
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	25 kDa
Uniprot ID	Q00059
Dilution Range	WB: 1:500-1:3000; IP: 1:500-1:1000; IHC: 1:50- 1:500
Expiration Date	12 months from date of receipt.