

Product Datasheet

TGM2 antibody (orb631858)

Description	Rabbit polyclonal antibody to TGM2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	transglutaminase 2(C polypeptide, protein-glutamine-gamma-glutamyltransferase)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	77 kDa, 82 kDa
Uniprot ID	P21980
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.