

Product Datasheet

TGN46 antibody (orb631861)

Description	Rabbit polyclonal antibody to TGN46
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	trans-golgi network protein 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	90-100 kDa
Uniprot ID	O43493
Dilution Range	WB: 1:1000-1:3000; IP: 1:200-1:1000; IF: 1:10- 1:100
Expiration Date	12 months from date of receipt.