

Product Datasheet

THG1L antibody (orb631872)

Description	Rabbit polyclonal antibody to THG1L
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	tRNA-histidine guanylyltransferase 1-like(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q9NWX6
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.