

Product Datasheet

CD71 antibody (orb632107)

Description	Rabbit polyclonal antibody to CD71
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	transferrin receptor(p90, CD71)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	77 kDa
Uniprot ID	P02786
Dilution Range	WB: 1:500-1:5000, IP: 1:200-1:2000, IHC-P: 1:100-1:400, IF: 1:50-1:500
Expiration Date	12 months from date of receipt.