

Product Datasheet

TrkB antibody (orb632164)

Description	Rabbit polyclonal antibody to TrkB
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	neurotrophic tyrosine kinase, receptor, type 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	43 kDa, Observed 55 kDa
Uniprot ID	Q16620
Dilution Range	WB: 1:500-1:5000, IP: 1:200-1:1000, IHC-P: 1:50- 1:500
Expiration Date	12 months from date of receipt.