

Product Datasheet

TRPC4AP antibody (orb632178)

Description	Rabbit polyclonal antibody to TRPC4AP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	transient receptor potential cation channel, subfamily C, member 4 associated protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	112 kDa
Uniprot ID	Q8TEL6
Dilution Range	WB: 1:200-1:1000, IP: 1:200-1:1000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.