

Product Datasheet

TRPV5 antibody (orb632190)

Description	Rabbit polyclonal antibody to TRPV5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	transient receptor potential cation channel, subfamily V, member 5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	98 kDa
Uniprot ID	Q9NQA5
Dilution Range	IF: 1:10-1:100
Expiration Date	12 months from date of receipt.