

Product Datasheet

TTC32 antibody (orb632236)

Description	Rabbit polyclonal antibody to TTC32
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	tetratricopeptide repeat domain 32
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	57 kDa
Uniprot ID	Q5I0X7
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:20-1:200, IF: 1:10-1:100, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.