

Product Datasheet

TTC7A antibody (orb632240)

Description	Rabbit polyclonal antibody to TTC7A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	tetratricopeptide repeat domain 7A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	45 kDa
Uniprot ID	Q9ULT0
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:20-1:200, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.