

Product Datasheet

TTLL10 antibody (orb632247)

Description	Rabbit polyclonal antibody to TTLL10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	tubulin tyrosine ligase-like family, member 10
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	43 kDa
Uniprot ID	Q6ZVT0
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.