

Product Datasheet

UBE2Q2 antibody (orb632315)

Description	Rabbit polyclonal antibody to UBE2Q2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	ubiquitin-conjugating enzyme E2Q family member 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	141 kDa, Observed 141 kDa, 200 kDa
Uniprot ID	Q8WVN8
Dilution Range	WB: 1:500-1:5000, ;IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.