

Product Datasheet

UBR2 antibody (orb632338)

Description	Rabbit polyclonal antibody to UBR2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	ubiquitin protein ligase E3 component n-recognin 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	53 kDa, Observed 70 kDa
Uniprot ID	Q8IWW8
Dilution Range	WB: 1:200-1:2000, IF: 1:10-1:100, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.