

Product Datasheet

PGP9.5 antibody (orb632345)

Description	Mouse monoclonal antibody to PGP9.5
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	ubiquitin carboxyl-terminal esterase L1(ubiquitin thiolesterase)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	57 kDa
Uniprot ID	P09936
Dilution Range	WB: 1:5000-1:50000, IF: 1:20-1:200, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.