

Product Datasheet

UCHL5IP antibody (orb632349)

Description	Rabbit polyclonal antibody to UCHL5IP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	UCHL5 interacting protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	37 kDa
Uniprot ID	Q99871
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:20-1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.