

Product Datasheet

ULK1 antibody (orb632378)

Description	Rabbit polyclonal antibody to ULK1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	unc-51-like kinase 1(C. elegans)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150 kDa
Uniprot ID	O75385
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.