

Product Datasheet

VAPA antibody (orb632484)

Description	Rabbit polyclonal antibody to VAPA
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	VAMP(vesicle-associated membrane protein)-associated protein A, 33kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	85 kDa
Uniprot ID	Q9P0L0
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:50-1:200;IF: 1:10-1:100, IP: 1:500-1:5000
Expiration Date	12 months from date of receipt.