

Product Datasheet

VGLUT1 antibody (orb632512)

Description	Rabbit polyclonal antibody to VGLUT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	solute carrier family 17(sodium-dependent inorganic phosphate cotransporter), member 7
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	31 kDa
Uniprot ID	Q9P2U7
Dilution Range	IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.