

Product Datasheet

VGLUT2 antibody (orb632513)

Description	Rabbit polyclonal antibody to VGLUT2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	solute carrier family 17(sodium-dependent inorganic phosphate cotransporter), member 6
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	62 kDa
Uniprot ID	Q9P2U8
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.