

Product Datasheet

Visfatin 1 antibody (orb632522)

Description	Rabbit polyclonal antibody to Visfatin 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	nicotinamide phosphoribosyltransferase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	52 kDa
Uniprot ID	P43490
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:20-1:200, IP: 1:500-1:2000
Expiration Date	12 months from date of receipt.