

Product Datasheet

Vitronectin antibody (orb632524)

Description	Rabbit polyclonal antibody to Vitronectin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	vitronectin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	53 kDa
Uniprot ID	P04004
Dilution Range	WB: 1:1000-1:10000, IF: 1:50-1:500, IHC-P: 1:100-1:400
Expiration Date	12 months from date of receipt.