

Product Datasheet

VMA21 antibody (orb632526)

Description	Rabbit polyclonal antibody to VMA21
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	VMA21 vacuolar H ⁺ -ATPase homolog(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	96 kDa, observed 160 kDa, 60 kDa
Uniprot ID	Q3ZAA7
Dilution Range	WB: 1:500-1:5000, IP: 1:500-1:5000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.