

Product Datasheet

WDR55 antibody (orb632599)

Description	Rabbit polyclonal antibody to WDR55
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	WD repeat domain 55
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	36 kDa, observed 42 kDa, 36 kDa
Uniprot ID	Q9H6Y2
Dilution Range	WB: 1:200-1:1000, IF: 1:10-1:100, IHC-P: 1:20-1:200, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.