



Product Datasheet

Endothelial lipase antibody (orb6326)


Description	Endothelial lipase antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Endothelial lipase (301-400/500 aa)
Target	Endothelial lipase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q9Y5X9
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat test...



Immunohistochemical staining of Rat lung...



Immunohistochemical staining of rat test...