

Product Datasheet

ZCCHC8 antibody (orb632705)

Description	Rabbit polyclonal antibody to ZCCHC8
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	zinc finger, CCHC domain containing 8
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	79 kDa, Observed 100 kDa
Uniprot ID	Q6NZY4
Dilution Range	WB: 1:500-1:5000, IF: 1:50-1:500
Expiration Date	12 months from date of receipt.