

Product Datasheet

ZIM3 antibody (orb632728)

Description	Rabbit polyclonal antibody to ZIM3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	zinc finger, imprinted 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	Q96PE6
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.