

Product Datasheet

ZNF311 antibody (orb632769)

Description	Rabbit polyclonal antibody to ZNF311
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	zinc finger protein 311
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	76 kDa, observed 70 kDa and 60 kDa
Uniprot ID	Q5JNZ3
Dilution Range	WB: 1:1000-1:10000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.