

Product Datasheet

ZNF839 antibody (orb632827)

Description	Rabbit polyclonal antibody to ZNF839
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	zinc finger protein 839
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	87 kDa
Uniprot ID	A8K0R7
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:20-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.