

Product Datasheet

ZO-1 antibody (orb632831)

Description	Rabbit polyclonal antibody to ZO-1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	tight junction protein 1(zona occludens 1)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	230kd
Uniprot ID	Q07157
Dilution Range	IF: 1:20-1:200
Expiration Date	12 months from date of receipt.