

Product Datasheet

ZO-2 antibody (orb632833)

Description	Rabbit polyclonal antibody to ZO-2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	tight junction protein 2(zona occludens 2)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	134 kDa, observed 150 kDa
Uniprot ID	Q9UDY2
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.