

Product Datasheet

Huntingtin antibody (orb632921)

Description	Rabbit polyclonal antibody to Huntingtin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Huntingtin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	348 kDa
Uniprot ID	P42858
Dilution Range	IHC-P: 1:50-1:500, IF: 1:50-1:500
Expiration Date	12 months from date of receipt.