

Product Datasheet

IgG control antibody (orb632959)

Description	Goat polyclonal antibody to IgG control
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	22 kDa
Dilution Range	WB: 1:2000-1:20000
Expiration Date	12 months from date of receipt.