

Product Datasheet

TUBB3 specific antibody (orb632967)

Description	Mouse monoclonal antibody to TUBB3 specific
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	AS160
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	18 kDa
Uniprot ID	Q13509
Dilution Range	WB: 1:5000-1:50000, IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.