

Product Datasheet

ATPB antibody (orb633020)

Description	Mouse monoclonal antibody to ATPB
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	ATP synthase, H ⁺ transporting, mitochondrial F1 complex, beta polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	56 kDa
Uniprot ID	P06576
Dilution Range	WB: 1:5000-1:30000, IHC-P: 1:200-1:1000
Expiration Date	12 months from date of receipt.