

Product Datasheet

VEGF Receptor 2-Y1175 antibody (orb633034)

Description	Rabbit polyclonal antibody to VEGF Receptor 2-Y1175
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1175 of human VEGF Receptor 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	250 kDa
Uniprot ID	P35968
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.