

Product Datasheet

YY1 antibody (orb633094)

Description	Rabbit polyclonal antibody to YY1
Species/Host	Rabbit
Reactivity	Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	YY1 transcription factor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa, 75 kDa
Uniprot ID	P25490
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.