

Product Datasheet

UTRO antibody (orb633151)

Description	Rabbit polyclonal antibody to UTRO
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Utrophin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P46939
Dilution Range	IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.