

Product Datasheet

MGMT antibody (orb633165)

Description	Rabbit polyclonal antibody to MGMT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	O-6-methylguanine-DNA methyltransferase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	22 kDa
Uniprot ID	P16455
Dilution Range	IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.