

Product Datasheet

CD100 antibody (orb633179)

Description	Mouse monoclonal antibody to CD100
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	sema domain, immunoglobulin domain(Ig), transmembrane domain(TM) and short cytoplasmic domain,(semaphorin) 4D
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	150 kDa
Uniprot ID	Q92854
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.