

Product Datasheet

6X His tag monoclonal antibody antibody (HRP) (orb633198)

Description	Mouse monoclonal antibody against 6X His tag monoclonal antibody conjugated to HRP
Species/Host	Mouse
Conjugation	HRP
Tested Applications	IP, WB
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:4000-1:10000 for Western Blotting with chromogenic substrates after immunoprecipitation
Expiration Date	12 months from date of receipt.