

Product Datasheet

Rabbit Anti-Monkey IgG antibody (HRP) (orb633271)

Description	Rabbit polyclonal antibody against Rabbit Anti-Monkey IgG conjugated to HRP
Conjugation	HRP
Tested Applications	ELISA, WB
Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Concentration	3mg/ml BSA (IgG-Free, Protease-Free)
Storage	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.