

Product Datasheet

Affinipure Rabbit Anti equine Ig (G, M, A) antibody (FITC) (orb633273)

Description	Rabbit polyclonal antibody against Affinipure Rabbit Anti equine Ig (G, M, A) conjugated to FITC
Conjugation	FITC
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	A minimum titer of 1:64 was determined by Fluorescent Dot Immunobinding Assay (FDIBA) on a 5 µg dot of equine Ig (G, M, A) bound to nitrocellulose. 1:100-1:200 for IF
Expiration Date	12 months from date of receipt.